




# Products & components



## Ceiling hoist system

### **GH3 Lifting modules** **4**

GH3	.5
GH3+	.5
GH3/GH3+ with drive motor	.5
GH3 Twin	.5
Optional modules for GH3+	.6
– CLM module – standard	
– CLM module – online	
– Scale module	
– Scale Class III module	
– Service module	
– Trainer module	
Technical specifications	14

### **GH1 Lifting modules** **8**

GH1	.9
GH1 Q	.9
GH1 F	10
Technical specifications	14

### **GHZ Lifting module** **12**

GHZ	13
Technical specifications	14

### **Rails** **15**

Rails	15
Wall rail	18
Components for rail systems	19
– Brackets	
– Trolleys	
– Components	
Free-standing rail systems	29

## Mobile lifters

GL5.2 Mobile lifter	30
GLS5.2 Active lifter	31
Variants – GL5.2 and GLS5.2	32
GTP1 – Transfer platform	33

## Slings and lifting accessories

ABC Slings – Active	36
ABC Slings – Basic	38
ABC Slings – Custom	40
ABC Slings – Kids	42
Bariatric slings	44
Disposable slings	47
Sizes and fabrics	52
Turning and handling options	56
Lifting hanger options	59
Horizontal lifters	61
Accessories	64



Transfer



Treatment



Training

*Guldmann is the expert at providing safe, comfortable and efficient solutions to improve the working environment in the healthcare sector where people need to be lifted, moved and positioned as well as mobilized and rehabilitated.*

*Our prime objective is to make sure that carers get more time to provide all-important personal care. We achieve this by enabling them to make all lifting and moving operations as streamlined and efficient as possible at every single stage.*

*And we do this in ways that meet the exact needs of all concerned. In this specialised field, experience really counts. Guldmann has perfected a comprehensive, modular range of lifting and moving solutions, drawing on more than a quarter of a century of experience in supplying assistive technology products designed to help people live their lives to the full.*

**The Guldmann approach**

*We develop and manufacture a wide range of products that take both the carer's working environment and the user's personal comfort and security into consideration. We help everyone involved to choose the solutions, products and modules that best meet their needs.*


*Our products are designed to be combined and used in a wide range of different ways, enabling us to provide you with a complete solution designed to match your current need as well as being well-aligned with the functionality requirements of tomorrow.*

*We provide analysis tools and expert advice to ensure a well-informed decision. We work with architects, consulting engineers and building contractors – and all the many other professions involved – to ensure that the solutions decided on can be integrated efficiently into both new and existing facilities.*

*We also help customers by ensuring these solutions are correctly installed. We run training programmes that help carers decide on the best solutions for their special needs, and train them in using the system in the safest, most effective way.*

*We are focused on creating more “time to care” in connection with lifting and moving requirements across the full spectrum of healthcare environments. Our services include support, education, installation and maintenance.*

**| Time to care |**



Ceiling hoist system


GH3 Lifting modules



## GH3 Lifting modules



red dot design award  
honourable mention 2009



GH3  
hand control

GH3+  
hand control

1

### GH3

Lifting module. User-friendly and extremely powerful and fast lifting module for permanent installation. Lifting capacity 200 or 250 kg

2



### GH3+

Lifting module similar to GH3, but with options like integrated scale, CareLift Management and Service monitoring system and even faster lifting speed. Lifting capacity 200, 250, 275, 300, 350, 375 or 400 kg.

3



### GH3/GH3+ with drive motor

For horizontal movement of the lifting module in the rail system the GH3 and GH3+ are available with drive motor. They can be operated with standard hand control (with cord) or with optional infrared (IR) remote control.

4



### GH3 Twin

Heavy duty lifting module with double lifting strap. Lifting capacity 250, 375 or 500 kg.

See table on page 14



5

**CareLift Management  
CLM module – standard**

Provides reliable data about the number and types of lift over a given period. This data is available as a downloadable spread for use in standard types of spreadsheet software.

6

**CareLift Management  
CLM module – online**


Provides instant overview and will keep track of when and how often any ceiling hoist is used. The CLM data from the hoist units are presented in the CLM Viewer.

7

**Service module**

Provides information about when the next scheduled safety inspection/service is due. This module also works with the CLM module online.





8

**Scale module**

Integrated scale module used for weighing of a person when lifted in a sling. The weight of the lifted person is shown on the hand control display.




9

**Scale Class III module**

The integrated Class III scale is used for weighing of patients in connection with e.g. hospitalisation where accurate monitoring of a patient's weight is necessary. The weight is displayed on the hand control.


This is the first Class III approved scale on the market that is integrated into a ceiling hoist fulfilling the conditions for non-automatic weighing instrument in accordance with EU Directive 2014/31/EU.



10

**Trainer module**

Provides rehabilitation with dynamic weight relief up to 100 kg in steps of 1 kg. Typical rehab exercises with dynamic weight relief in the Trainer Module could be: Gait training, balance, sit-to-stand, squat exercises, steps on/over obstacles etc.



Ceiling hoist system

GH1 Lifting modules





GH1 Lifting modules



11

**GH1**

Lifting module. User-friendly, powerful and fast lifting module for permanent installation and for installation in free-standing rail system. Lifting capacity 175, 205 or 255 kg. See table on page 14

12



**GH1 Q**

GH1 Q lifting module with a quick release system which makes it easy to click on and click off the lifting module from the rails for relocating. Lifting capacity 175, 205 or 255 kg. *GH1Q must not be used in systems with Combi-lock or switch track. See table on page 14*

13



**GH1 covers**

The GH1 lifting module is delivered standard with a choice of two options.

14



GH1 hand control

GH1 docking station

Ceiling hoist system

GH1 F Lifting module



**GH1 F Lifting module**



15

**GH1 F**

Flexible lifting module that can be moved from one rail system to another. It installs and removes itself automatically from the rail system. Lifting capacity 175, 205 or 255 kg. See table on page 14

16



**Travelling trolley for GH1 F**  
.....551165

17



**Transport trolley for GH1 F**  
GH1 F transport and storage trolley.  
.....554273

18



*GH1 F hand control*

*GH1 F docking station*

Ceiling hoist system

GHZ Lifting module



GHZ Lifting module



19



**GHZ**

A discreet and almost invisible ceiling hoist system. The GHZ lifting module runs inside the traverse rail. This minimizes the built-in dimensions and adds to the lifting height which is often a crucial factor in rooms with lower ceiling heights. Mounted in the wall

rails the flexible GHZ traverse rail will automatically compensate for non-aligned (oblique) walls with up to 80 mm. Maximum lifting capacity: 255 kg  
*See table on page 14*

.....556342

20




*GHZ hand control*

*GHZ docking station*

Technical specifications

GH Lifting modules, configurations						Options												
Guldmann hoist type	Product line	Load in kg	Number of lifting straps	Number of lifting motors	Number of horizontal drive motors <sup>1)</sup>	IR remote <sup>2)</sup>	Service module	CLM module	Scale module	WiFi module	Class III scale	Trainer module <sup>2)</sup>	Coating	Excl. Turbospeed	UL	Charg. HC <sup>2)</sup>		
GH3	(x)	xxx	x	x	x	Configuration code (xxxxx)												
GH3		200	1	1	0 - 1 - 2	•										•		
		250	1	1	0 - 1 - 2	•										•		
	+	200	1	1	0 - 1 - 2	•	•	•	•	•	•	•	•	•	•	•	•	•
		250	1	1	0 - 1 - 2	•	•	•	•	•	•	•	•	•	•	•	•	•
		275	1	1	0 - 1 - 2	•	•	•	•	•	•	•	•	•	•	•	•	•
		300	1	2	0 - 2	•	•	•	•	•	•	•	•	•	•	•	•	•
		350	1	2	0 - 2	•	•	•	•	•	•	•	•	•	•	•	•	•
		375	1	2	0 - 2	•	•	•	•	•	•	•	•	•	•	•	•	•
		400	1	2	0	•	•	•	•	•	•	•	•	•	•	•	•	•
	Twin	250	2	2	0	•	•	•	•	•	•			•	•	•		
		375	2	2	0	•	•	•	•	•	•			•	•	•		
		500	2	2	0	•	•	•	•	•	•			•	•	•		
GH1		175	1	1	0									•		•		
		205	1	1	0									•		•		
		255	1	1	0									•		•		
	F	175	1	1	0										•		•	
		205	1	1	0										•		•	
		255	1	1	0										•		•	
	Q	175	1	1	0										•		•	
		205	1	1	0										•		•	
		255	1	1	0										•		•	
GHZ		255	1	1	0									•		•		

1) Drive motors are not compatible with UL  
 2) Not compatible with UL



## Rails



1

Rails A/B/C								
Rail A	1 m	2 m	3 m	4 m	5 m	6 m	7 m	8 m
w/ charging liner	550351	550352	550353	550354	550355	550356		
Cut-out for trolley	552535	552536	552537	552538	552539	552540		
w/o charging liner	550670	550671	550672	550673	550674	550675		
Cut-out for trolley	557882	557883	557884	557885	557886	557887		
Rail B	1 m	2 m	3 m	4 m	5 m	6 m	7 m	8 m
w/ charging liner			550359	550360	550361	550362	550363	550364
Cut-out for Combi-lock – short			555620	555621	555622	555623	555624	555625
Cut-out for Combi-lock – long			555614	555615	555616	555617	555618	555619
Cut-out for trolley			552541	552542	552543	552544	552545	552546
w/o charging liner			550677	550678	550679	550680	550681	550682
Cut-out for trolley			557888	557889	557890	557891	557892	557893
Rail C	1 m	2 m	3 m	4 m	5 m	6 m	7 m	8 m
w/ charging liner					550365	550366	550367	550368
Cut-out for Combi-lock – short					555630	555631	555632	555633
Cut-out for Combi-lock – long					555626	555627	555628	555629
Cut-out for trolley					552547	552548	552549	552550
w/o charging liner					550683	550684	550685	550686
Cut-out for trolley					557894	557895	557896	557897

The ABC rails are also available in white (RAL 9016). If white is requested add -WH to the item no. (xxxxx-WH).

### Special rail cuts

Cut-out for Combi-lock  
See table above

Cut-out for trolley  
See table above



Rails



2

Rails A+/B+/C+								
Rail A+	1 m	2 m	3 m	4 m	5 m	6 m	7 m	8 m
w/ charging liner	562361	562362	562363	562364	562365	562366		
Cut-out for trolley	562377	562378	562379	562380	562381	562382		
w/o charging liner	562601	562602	562603	562604	562605	562606		
Cut-out for trolley	562617	562618	562619	562620	562621	562622		
Rail B+	1 m	2 m	3 m	4 m	5 m	6 m	7 m	8 m
w/ charging liner			562367	562368	562369	562370	562371	562372
Cut-out for trolley			562383	562384	562385	562386	562387	562388
w/o charging liner			562607	562608	562609	562610	562611	562612
Cut-out for trolley			562623	562624	562625	562626	562627	562628
Rail C+	1 m	2 m	3 m	4 m	5 m	6 m	7 m	8 m
w/ charging liner					562373	562374	562375	562376
Cut-out for trolley					562389	562390	562391	562392
w/o charging liner					562613	562614	562615	562616
Cut-out for trolley					562629	562630	562631	562632

The ABC rails are also available in white (RAL 9016). If white is requested add -WH to the item no. (xxxxx-WH).

3

Cut-out for trolleys 250 mm

	Grey	White
Rail A	558926	558926-WH
Rail B	558927	558927-WH
Rail A+	562337	562337-WH
Rail B+	562338	562338-WH

4


Cut-out for trolleys 620 mm

	Grey	White
Rail A	558928	558928-WH
Rail B	558929	558929-WH
Rail A+	562339	562339-WH
Rail B+	562340	562340-WH





Rails



5

**Rail curves (rail A)**

	Grey	White
30° .....	550015 .....	550015-WH
45° .....	550016 .....	550016-WH
60° .....	550017 .....	550017-WH
90° .....	550018 .....	550018-WH
90° large .....	550025 .....	550025-WH



6

**Rail curves (rail A+)**

	Grey	White
30° .....	562409 .....	562409-WH
45° .....	562410 .....	562410-WH
60° .....	562411 .....	562411-WH
90° .....	562412 .....	562412-WH
90° large .....	562413 .....	562413-WH

Rails



7

8




**Wall rail**

	Rail	White cover
3 m. ....	560813. ....	556013
4 m. ....	560814. ....	556014
5 m. ....	560815. ....	556015
7 m. ....	560817. ....	556017

**Wall rail with charging liner**

	Rail	White cover
3 m. ....	560803. ....	556013
4 m. ....	560804. ....	556014
5 m. ....	560805. ....	556015
7 m. ....	560807. ....	556017



9



**GHZ traverse rail**

	White
3 m. ....	557703
4 m. ....	557704
5 m. ....	557705

Components for rail systems



10

**Combi-lock, 375 kg**

..... 555237

11



**Turntable, electric**

	Grey	White
255 kg .....	550663 .....	550663-WH .....
- with coating .....	551430 .....	551430-WH .....
375 kg .....	552882 .....	552882-WH .....
- with coating .....	552984 .....	552984-WH .....

12



**Switch track, electric, 250 kg**

	Grey	White
60° right .....	553991 .....	553991-WH .....
- with coating .....	551058 .....	551058-WH .....
60° left .....	553995 .....	553995-WH .....
- with coating .....	551059 .....	551059-WH .....
90° right .....	553925 .....	553925-WH .....
- with coating .....	554027 .....	554027-WH .....
90° left .....	553926 .....	553926-WH .....
- with coating .....	554028 .....	554028-WH .....



13

**Rail crossing module**

Flexible division for set-ups with mobile walls  
and folding doors .....

554528

Components for rail systems

14



**Wall bracket end**  
*end-mounted rail*  
 Grey ..... 550451  
 White ..... 550471

15



**Wall bracket end**  
*side-mounted rail*  
 Grey ..... 550451  
 White ..... 550471

16



**Wall bracket side**  
*end-mounted rail*  
 Grey ..... 550450  
 White ..... 550472

17



**Wall bracket side**  
*side-mounted rail*  
 Grey ..... 550450  
 White ..... 550472

18



**Rail joint bracket**  
 ..... 12448

19



**Ceiling bracket, standard**  
**for straight rail and curve**  
 Grey ..... 550402  
 White ..... 550433

20



**Ceiling bracket, standard**  
**for rail joints**  
 Grey ..... 550402  
 White ..... 550433

21



**Assembly bracket**  
 ..... 12336

22



**Ceiling bracket adjustable 5 mm**  
**for Combi-lock**  
 Grey ..... 555330  
 White ..... 555567

23



**Lowered ceiling bracket standard**  
**for straight rail and curve**  
 Grey ..... 550403  
 White ..... 550434

24



**Plasterboard bracket**  
**max 255 kg**  
 White ..... 550452

Components for rail systems

25



**Ceiling bracket, straight lowered for straight rail and curve**

Length	Grey	White
250-450 mm	562015	562015-WH
451-650 mm	562016	562016-WH
651-850 mm	562017	562017-WH
851-1000 mm	562018	562018-WH
1001-1200 mm	562019	562019-WH
1201-1400 mm	562020	562020-WH
1401-1600 mm	562021	562021-WH
1601-1800 mm	562022	562022-WH
1801-2000 mm	562023	562023-WH

26



**Ceiling bracket, straight lowered for Combi-lock**

Length	Grey	White
316-520 mm	562060	562060-WH
521-720 mm	562061	562061-WH
721-920 mm	562062	562062-WH
921-1125 mm	562063	562063-WH
1071-1270 mm	562064	562064-WH
1271-1470 mm	562065	562065-WH
1471-1670 mm	562066	562066-WH
1671-1870 mm	562067	562067-WH
1871-2070 mm	562028	562028-WH

27



**Upright support brackets**

*White as standard, RAL 9016*

without foot	124841
with low foot (150 mm)	124801
with high foot (250 mm)	124811
Guideway for upright support	1248031

*Grey is an option, RAL 9006*

without foot	124844
with low foot (150 mm)	124804
with high foot (250 mm)	124814
Guideway for upright support	1248034

28



**Ceiling bracket, lowered rod for straight rail and curve**

Length	Grey	White
250 mm	550498	550499
500 mm	550436	550410
1000 mm	550437	550408
2000 mm	550438	550409

29



**Stabilizing bracket**

Length	Grey	White
500 mm	550417	550440
1000 mm	550418	550441
2000 mm	550419	550442

30



**Ceiling bracket slope lowered**

Length	Grey	White
250 mm	550412	550476
500 mm	550413	550477
1000 mm	550414	550474
2000 mm	550415	550475

Components for rail systems

31



**Traverse trolley**

For room covering systems mounted below the ceiling. Max. SWL 375 kg.

Grey . . . . .550076  
White . . . . .550491

32



**Traverse trolley Lowered 25 mm**

For room covering systems where parallel rails are integrated/recessed into the ceiling. Max SWL 375 kg.

Grey . . . . .554189  
White . . . . .551327

33



**Traverse trolley cross Combi-lock**

Must be used in room covering systems on the parallel rail closest to the Combi-lock or in systems with overhang of traverse rail. Max. SWL 375 kg.

Grey . . . . .550779  
White . . . . .555205

34



**Traverse trolley cross Combi-lock, lowered 25 mm**

Must be used in room covering systems on the parallel rail that is integrated/recessed into the ceiling and closest to the Combi-lock or in systems with overhang of traverse rail. Max. SWL 375 kg.

Grey . . . . .554526  
White . . . . .558664

35



**Traverse trolley**

For room covering systems mounted below the ceiling. Max. SWL 500 kg.

Grey . . . . .550078  
White . . . . .550493

36



**Traversløbekat Heavy Duty, forlænget 25 mm**

For room covering systems where parallel rails are integrated/recessed into the ceiling. Max. SWL 500 kg.

Grey . . . . .554546  
White . . . . .554548

37



**Traverse trolley Heavy duty cross**

For room covering systems with overhang of traverse rail. Max. SWL 500 kg.

Grey . . . . .559548  
White . . . . .559549

38



**Traverse trolley long**

For room covering systems where the traverse rail is longer than 5 meter. Max SWL 400 kg.

Grey . . . . .556106  
White . . . . .555802

39



**Traverse trolley long, lowered 25 mm**

For room covering systems where the parallel rails are integrated/recessed into the ceiling. Max SWL 400 kg.

Grey . . . . .558102  
White . . . . .560080

40



**Traverse trolley for wall rails B traverse rail**

Can only be used in combination with wall rail and B rail (traverse).

White . . . . .560790  
Kit for C-traverse rail . . . . .560720

Components for rail systems



**41 GHZ traverse trolley**  
Can only be used in combination with wall rails and GHZ traverse rail.  
White .....556290



**42 Adjustable traverse trolley lowered 25 mm**  
Max SWL 375 kg  
Grey .....561091  
White ..... 561091-WH



**43 Traverse trolley short, lowered 25 mm**  
To be used when the hoist has to be placed/stored in a cupboard. Max SWL 375 kg  
*Set of two (left and right)*  
Grey .....561960  
White ..... 561960-WH



**44 Drive motor for traverse rail**  
for traverse B and C rails. ....562425



**45 Trolley assembly GH2 F and GH1 F – for curtain**  
Trolley that can pass through a rail gab (curtain rail for privacy).  
..... 12776-8



**46 Trolley GH3 – for curtain**  
Trolley that can pass through a rail gab (curtain rail for privacy).  
.....550581



**47 Trolley GH3+ with loadcell for curtain**  
Trolley that can pass through a rail gab (curtain rail for privacy).  
.....550582



**48 Trolley assembly ext. 25 mm GH3**  
.....554081



**49 Trolley assembly ext. 50 mm GH3**  
.....552503



**50 Pick-up assembly traverse trolley**  
For A, B, C and A+, B+, C+ rails  
.....550094



**51 Pick-up assembly for Wall rail**  
.....560770



**52 Kit for C-rail – Wall rail**  
.....560720

Components for rail systems

53



**Trolley assembly ext. 25 mm for GH1**  
.....554082

54



**Trolley drive assembly ext. 25 mm  
1 drive**  
.....554263

55



**Travelling trolley for GH1 Q**  
.....552730

56



**Rail joint tool**  
.....561164

57



**Connector plate – wall rail**  
.....561446

58



**Profile for ceiling, end**  
.....559229

59

**Spacer lowering traverse rail**



	Grey	White
15 mm	.....551077	..... 551077-WH
20 mm	.....551078	..... 551078-WH
30 mm	.....551079	..... 551079-WH
35 mm	.....551096	..... 551096-WH
40 mm	.....551080	..... 551080-WH
50 mm	.....551081	..... 551081-WH
60 mm	.....551082	..... 551082-WH
70 mm	.....551083	..... 551083-WH
80 mm	.....551084	..... 551084-WH
90 mm	.....551085	..... 551085-WH
100 mm	.....551086	..... 551086-WH

*Each package includes bolts and spacers*





Components for rail systems



60

Side cover		1 m	2 m	3 m	4 m	5 m	6 m	7 m	8 m	100 m
<b>Rail A side cover</b>	Grey	550230	550231	550232	550233	550234	550235			550763
	White	550238	550239	550240	550241	550242	550243			550766
<b>Rail B side cover</b>	Grey			550246	550247	550248	550249	550250	550251	550764
	White			550252	550253	550254	550255	550256	550257	550767
<b>Rail C side cover</b>	Grey					550258	550259	550260	550261	550765
	White					550262	550263	550264	550265	550768

61



**Cover for wall rail**

	White
3 m. ....	556013
4 m. ....	556014
5 m. ....	556015
7 m. ....	556017
100 m. ....	556078

62



**End cover**

	Grey	White
Rail A .....	550447	550463
Rail B .....	550448	550464
Rail C .....	550449	550465

63



**End cover, white – wall rail**

.....	560776
-------	--------

Components for rail systems

64



**End stop cover – wall rail**  
.....560773

65



**Top cover**  
50 m. ....561411  
100 m. ....561410


66



**Cover for trolley**  
Grey. ....550570  
White .....550584



67



**Sealing Kit**  
1 set. ....559720

68



**Accessory kit – Wall rail**  
2 sets. ....560780

69



**Charging liner connector**  
.....551432

70



**Charging liner**  
6 m rail. ....550425  
10 m rail. ....550427

71

**Charging liner – wall rail**  
4 m .....561447  
10 m. ....561448

Components for rail systems

72



**Hand control for extension cable**  
 GH3+ .....552300  
 GH3+ drive .....552307  
 GH3 .....552302  
 GH3 drive .....552301  
 GH1 .....552619  
 GH1 F .....555903

73



**Extension cable for hand control**  
 GH3+ .....552303  
 GH3 .....552304  
 GH1 .....552620  
 GH1 F .....555902

74



**Remote control IR**  
 for switch track/turntable/GH3. . . . .935034  
 for Positioning lock . . . . .935039

**Hand control with cable**  
 for switch track/turntable. . . . .9350151

75



**Wall switch – flush mounted**  
 for switch track/turntable. . . . .9350471

**Wall switch – surface mounted**  
 for switch track/turntable. . . . .935049

76



**Power supply**  
 Transformer AC, Class I . . . . .550200  
 Power supply, GH1 . . . . .551244

77



**Cable main power AC – white**  
 EU .....554250  
 USA .....554251  
 UK .....554252  
 AU .....554253  
 CAN.....555878


Components for rail systems

78



**End stop**  
 Grey ..... 12330  
 White ..... 123301

79



**Set of brakes for activity rail**  
 ..... 121664

80



**Suspension eye for activity rail**  
 ..... 121662

81



**Positioning lock – single**  
*Must not be used in systems with Combi-lock*  
 ..... 562421

82



**Positioning lock – room covering**  
*Must not be used in systems with Combi-lock*  
 Positioning lock – room covering ..... 562423  
 Positioning lock – complete system (room covering and single) ..... 562424

83



**Weight cart, electrical**  
 Incl. weights, 550 kg ..... 290880

84



**Weight slice, 16.7 kg for manual weight cart**  
 ..... 558468

85



**Weight cart, manual**  
 Max 455 kg  
 incl. 2 weight slices, 49.5 kg ..... 290851  
 incl. 17 weight slices, 300 kg ..... 290850

Free-standing rail systems



86

**Free-standing room covering rail system  
– height and length adjustable**

Max capacity 255 kg

Height min 2040 mm / max 2840 mm<sup>x)</sup>

3 x 3 m (length min 3200 mm / max 4200 mm<sup>x)</sup> . . . . . 554729

4 x 4 m (length min 4200 mm / max 5200 mm<sup>x)</sup> . . . . . 557968

87



**Swing lift II**

Max lifting capacity 205 kg

Length 2.5 m

..... 17200

88



**Free-standing single rail system  
– height and length adjustable**

Max capacity 255 kg

Height min 2020 mm / max 2820 mm<sup>x)</sup>

2.0 m (length min 2220 mm / max 3390 mm<sup>x)</sup> . . . . . 556000

2.5 m (length min 2720 mm / max 3890 mm<sup>x)</sup> . . . . . 553000

3.0 m (length min 3220 mm / max 4390 mm<sup>x)</sup> . . . . . 556297

<sup>x)</sup> External measurements

GL5.2 Mobile lifter



1



**GL5.2**  
Mobile lifter with electric width adjustment.  
Max lifting capacity 155 or 205 kg.  
*See table on page 32*

GLS5.2 Active lifter



2



**GLS5.2**

Active lifter with electric width adjustment.  
Max. lifting capacity 155 or 205 kg  
*See table on page 32*

3



**GLS5.2 Standing sling**

Supports the upper back under the shoulders to lower back.

**Suitable for**

- persons who can stand up but not raise themselves.

S/M	.....	28336	.....	(belt M	... 281553)
L	.....	28337	.....	(belt XL	.. 281573)
XL	.....	28339	.....	(belt 2XL	. 281583)
2XL	.....	28340	.....	(belt 3XL	. 281593)

4



**GLS5.2 Hip support**

Supports the hip part of the body.

**Suitable for**

- persons who can stand up but have a tendency to let the hip slide back.

.....28338

Variants – GL5.2 and GLS5.2

Item no.	Item name	Castors 100 mm	Castors 80 mm	Manuel lowering
562308	GL5.2 205 480 0000	x		
562309	GL5.2 155 480 0000	x		
562310	GL5.2 205 380 0000		x	
562311	GL5.2 155 380 0000		x	
562312	GL5.2 205 460 0000	x		x
562313	GL5.2 155 460 0000	x		x
562314	GL5.2 205 360 0000		x	x
562315	GL5.2 155 360 0000		x	x
552200	GLS5.2 155 560 0000			x
556001	GLS5.2 205 560 0000			x



5

6



**Battery pack**  
 Easy exchange and external charging possible  
 .....561553

**External charger**  
 .....561556



GTP1 – Transfer platform



7

8



**GTP1**

The GTP1 is a manual, mobile working tool that makes it easier to move standing and semi-independent users. Max. weight capacity of 175 kg.

.....560700

**GTP1 Transfer belt**

Supports the upper body.

**Suitable for**

- persons who can stand up themselves but need extra support during transfer.

S .....	284141
M .....	284151
L .....	284161



# Guldmann slings – a host of options


*Using the appropriate sling is essential for efficient care, and to ensure comfort, dignity and safety for the patient without compromising the safety of the care provider.*

*You can find a sling for almost any care task in the extensive range of slings and lifting accessories, which can be used with both ceiling hoists and mobile lifters.*

Transfer



Positioning



Lifting



Treatment



Mobilisation  
and activation



Training



ABC Slings – Active

1



2



3



**Active Trainer**

Provides support around the chest and shoulder areas. Detachable leg straps give support around the hip area and prevent the user from sliding out of the sling.

**Suitable for**

- persons with reduced body balance, but must be able to bear weight on their legs
- walking training • supported standing and upright postures • early mobilization.

..... 2830xx

**Active Micro**

Supports the upper part of the body (from the pelvis to under the shoulders) around the chest and upper thigh. The chest straps tighten around the chest.

**Suitable for**

- lifting persons with head control, but with a reduced muscular strength of upper part of the body and above hip joint and thighbone
- for some persons with extension spasms
- in situations where the sling should be put on and taken off easily
- lifting to and from a seated position
- situations where a large opening at the bottom is required, such as access for toileting and hygiene.

..... 2840xx

**Active Micro Plus**

Supports the upper part of the body (from the pelvis to under the shoulders) around the chest and upper thigh. The chest straps tighten around the chest.

**Suitable for**

- lifting persons with head control, but with a reduced muscular strength of upper part of the body and above hip joint and thighbone
- for some persons with extension spasms
- in situations where the sling should be put on and taken off easily
- lifting to and from a seated position
- situations where a large opening at the bottom is required, such as access for toileting and hygiene.

..... 281xxx



## ABC Slings – Active



4

### Gait Trainer

Supports the upper part of the body and has leg cuffs for enhanced safety. Added comfort due to Spacer material.

### Suitable for

- persons who have sufficient leg strength to stand upright, but need support due to reduced balance
- walking training
- supported standing and upright postures
- balance training early mobilization.

The sling is among other things designed for use in conjunction with the Guldmann Trainer Module.

..... 2832xx



ABC Slings – Basic

5



6



7



**Basic Basic**

Supports the whole body excluding the head.

**Suitable for**

- lifting persons with head control and with or without body control
- lifting to and from a seated or half-seated position (bed and wheelchair etc)
- lifting and transfers from all patient furniture surfaces.

..... 270xxx

**Basic Low**

Supports the body up to shoulder height.

**Suitable for**

- lifting persons with head and body control
- lifting to and from a seated position, (chair and bed etc.)
- lifting and transfers from all patient furniture surfaces.

..... 271xxx

**Basic High**

Supports the whole body including the head.

**Suitable for**

- lifting persons with little control over their head and body
- lifting to and from a lying position (bed and floor etc.)
- lifting to and from a seated position
- lifting and transfers from all patient furniture surfaces and the floor.

..... 272xxx



ABC Slings – Basic

8



9



**Basic Shell**

Supports the whole body and the head, but allows the arms to be outside the sling.

**Suitable for**

- lifting persons with little control over their head and body
- lifting to and from a lying position (bed and floor etc.)
- lifting to and from a seated position – letting the sling stay in the wheelchair.

..... 2750xx

**Basic Comfort High**

Supports the whole body and the head. The adjustable head support can be folded down when user is left in sitting position on the sling. Adds extra support around the hip and thigh area. Can be used as divided legs or as a cross over position, this also provides good support for the amputee.

**Suitable for**

- lifting persons with little control over their body
- lifting to and from a lying position (bed and floor etc.)
- lifting to and from a seated position
- persons who need added support around the hip and thighs + in situations where it is important that the thigh bone remains in a neutral position.

..... 2770xx



◀ Watch instruction video

ABC Slings – Custom

10



11



12



**Sit-On II**

Support the thighbone and back. The user sits on the sling. The design of the sling prevents the thighbone from an invert rotation when lifting.

**Suitable for**

- persons who need to stay sitting on the sling for long periods due to high weight or customized wheelchair.

..... 2970xx

**Sit-On High II**

Supports the thighbone, back and head. The design prevents the thighbone from an invert rotation when lifting. The head support can be adjusted to suit the user's needs and can be folded away whilst sitting in a wheelchair.

**Suitable for**

- persons who need to stay sitting on the sling for long periods due to high weight or customized wheelchair.

..... 2980xx

**Custom Amputee**

Supports the whole body up to shoulder height and the pelvis region.

**Suitable for**

- lifting both single and double (above knee joint) amputees
- lifting double and high amputees
- lifting from a lying position. Head support should be used to support head.

..... 2900xx





ABC Slings – Custom

13



14



15



**Sit-On Comfort**

Supports the thighbone and back. The user sits on the sling. The design of the sling prevents the thighbone from an invert rotation when lifting. Added comfort due to the Spacer material.

**Suitable for**

- persons who remain on the sling for long periods due to high weight or customized wheelchair.

*Other features*

• material prevents heat generation between body and sling • sling is fitted on the person while lying, but is used when lifting to and from a seated position in bed or wheelchair.

..... 2930xx

**Sit-On Comfort High**

Supports the thighbone, back and head. The adjustable head support can be folded down when user is left in sitting position on the sling. The design of the sling prevents the thighbone from an invert rotation when lifting. Added comfort due to the Spacer material.

**Suitable for**

- persons who remain on the sling for long periods due to high weight or customized wheelchair.

*Other features*

• material prevents heat generation between body and sling • sling is fitted on the person while lying, and is used when lifting to and from a lying position.

..... 2940xx

**Sit-On Comfort High Hygiene**

Designed with an opening in the bottom of the sling. Supports the thighbone, back and head. The adjustable head support can be folded down. The design of the sling prevents the thighbone from an invert rotation when lifting. Added comfort due to the Spacer material.

**Suitable for**

- transfer from bed to toilet or toilet chair. It is not recommended to use the Sit-On Comfort High Hygiene sling as a traditional Sit-On sling in a wheelchair or similar device due to the risk of pressure on the sling's opening.

*Other features (see Sit-On Comfort High)*

..... 2941xx



ABC Slings – Kids

16



17



18



**Basic Basic – Kids**

Supports the whole body excluding the head. The inside belt offers an extra support of the body.

**Suitable for**

• lifting kids with head control and with or without body control • lifting to and from a seated or half-seated position (bed and wheelchair etc.).

..... 270xxx

**Basic High – Kids**

Supports the whole body including the head. The inside belt offers an extra support of the body.

**Suitable for**

• lifting kids with little control over their head and body • lifting to and from a lying position (bed, changing bench and floor etc.) • lifting to and from a seated position.

..... 272xxx

**Active Vest – Kids**

Supports the body and pelvis. The width can be adjusted.

**Suitable for**

• kids with reduced body balance • walking training • supports standing and upright postures • playing.

..... 2831xx



ABC Slings – Kids

19



20



21



**Sit-On Comfort – Kids**

Supports the thighbone and back. The user sits on the sling. The design of the sling prevents the thighbone from an invert rotation when lifting. Added comfort due to the Spacer material.

**Suitable for**

- kids who remain on the sling for long periods due to customized wheelchair, or need lifting regularly.

*Other features*

- material prevents heat generation between body and sling
- sling is fitted on the child while lying, but is used when lifting to and from a seated position in bed or wheelchair

..... 2930xx

**Sit-On Comfort High – Kids**

Supports the thighbone, back and head. The adjustable head support can be folded down when user is left in sitting position on the sling. The design of the sling prevents the thighbone from an invert rotation when lifting. Added comfort due to the Spacer material.

**Suitable for**

- kids who remain on the sling for long periods due to high weight or customized wheelchair.

*Other features*

- material prevents heat generation between body and sling
- sling is fitted on the child while lying, and is used when lifting to and from a lying position.

..... 2940xx

**Vest for Standing Shell – Kids**

Supports the head and back

**Suitable for**

- lifting the child from a lying to a standing position
- used in conjunction with a standing shell to assist into a standing frame or table.

..... 2835xx

Slings and lifting accessories

Bariatric slings



Bariatric slings

22



**Gait Trainer Bariatric, 500 kg**  
Supports the pelvis area and pannus.

**Suitable for**

- supporting the patient in a standing position, maximising comfort and support.
- One size . . . . .283100

23



**Basic High Bariatric, 500 kg**  
Supports the whole body including the head.


**Suitable for**

- lifting persons with little control over their head and body
  - lifting to and from a lying position (bed and floor etc.)
  - lifting to and from a seated position
  - lifting and transfers from all patient furniture surfaces and the floor.
- 3XL standard . . . . . 2720913  
4XL standard . . . . . 2720914  
5XL standard . . . . . 2720916



◀ Watch instruction video

Bariatric slings



24



**Leg sling**



Used when legs need lifting and support. Possible activities are dressing, apply therapeutic sock/hose, access to hygiene, support of leg to apply ROM devices.  
Supplied in a box of 10 slings. . . . . 28650

25



**Pannus Support**

Supports and lifts the pannus in a position that will ease the personal hygiene task and avoiding a manual lift of the pannus.

**Suitable for**

- bariatric persons with a large pannus
  - personal hygiene tasks in bed.
- ..... 28660



26



**Twin Turner Bariatric**

Makes it possible to roll a patient on both sides without removing the sling. Enables placing of sling under patients without having to roll or move them manually. Smaller surface area of material for this sling allows access to patient's skin for inspection, treatment and care.

..... 28760

27



**Repositioning sling Bariatric, 500 kg**

A sling for use as part of the linen make-up of the patient bed; leave positioned underneath patient, ready to use for care tasks.

**Suitable for**

- rolling to side lying
- repositioning patient up in bed
- sitting forward head and trunk
- supine lateral transfers to another surface such as beds + stretchers, full lift in supine and for weighing the patient, stripping bed linen, or other.

Supplied in a box of 5 slings. . . . . 284656



Disposable slings



*For efficient and hygienic lifting routines*

*Guldmann disposable slings are ideal for situations in which an extra-high hygiene standard is required due to risk of infection. The disposable slings are intended for single-patient use only and should be disposed of after discharge or after use when soaked or soiled.*

## Disposable slings

28



29



30



### Disposable High, 205 kg

Supports the whole body including the head. Same design as Guldmann Basic High sling.

#### Suitable for

- lifting persons without head and body control
- lifting to and from a lying position (bed, floor, etc)
- lifting to and from a seated position.

#### Material

Polypropylene

Supplied in a box of 10 slings. . . . . 2836xx

### Disposable High – Kids, 205 kg

Supports the whole body including the head. Same design as Guldmann Basic High – Kids sling.

#### Suitable for

- lifting kids without head and body control
- lifting to and from a lying position (bed, floor, etc)
- lifting to and from a seated position.

#### Material

Polypropylene

Supplied in a box of 10 slings. . . . . 2836xx

### Disposable High Bariatric, 500 kg

Supports the whole body including the head. Same design as Guldmann Basic High sling.

#### Suitable for

- lifting persons without head and body control
- lifting to and from a lying position (bed, floor, etc)
- lifting to and from a seated position.

#### Material

Polyester/cotton

Supplied in a box of 10 slings. . . . . 2836xx



## Disposable slings

31



32



33



### Disposable Comfort High, 255 kg

Supports the whole body. Adds extra support around the hip and thigh area. Can be used as divided legs or as a cross over position, this also provides good support for the amputee. Same design as Guldmann Basic Comfort High sling (see page 39).

#### Suitable for

- lifting persons with little control over their body
- lifting to and from a seated position
- persons who need added support around the hip and thighs + in situations where it is important that the thigh bone remains in a neutral position.

#### Material

Polyester  
Supplied in a box of 5 slings. . . . . 2770xx

### Disposable Micro Plus, 205 kg

Supports the upper part of the body (from the pelvis to under the shoulders) around the chest and upper thigh. The chest straps tighten around the chest.

#### Suitable for

- lifting persons with head control, but with a reduced muscular strength of upper part of the body and above hip joint and thighbone
- for some persons with extension spasms
- in situations where the sling should be put on and taken off easily
- lifting to and from a seated position – situations where a large opening at the bottom is required, such as access for toileting and hygiene.

#### Material

Polyester  
Supplied in a box of 5 slings. . . . . 2815xx

### Disposable Gait Trainer, 205 kg

Provides support around the chest and shoulder areas. Leg straps give support around the hip area and prevent the user from sliding out of the sling.

#### Suitable for

- persons with reduced body balance but able to carry weight on their legs
- walking training
- supported standing and upright postures
- early mobilization.

#### Material

Polyester  
Supplied in a box of 5 slings. . . . . 2835xx

Disposable slings

34



**Disposable  
Repositioning sling, 500 kg**

A sling for use as part of the linen make-up of the patient bed; leave positioned underneath patient, ready to use for care tasks. Same design as Guldmann Repositioning sling.

**Suitable for**

- rolling to side lying, repositioning patient up in bed
- sitting forward head and trunk
- supine lateral transfers to another surface such as beds + stretchers
- full lift in supine off of bed to zero bed scale
- change linen, or other.

**Material**

Polyester  
Standard (photo) . . . . . 284555  
Supplied in a box of 10 slings

35



**Disposable OR**

Supports the upper body when moving to and from the operating table in a horizontal position.

**Suitable for**

- the length and design of the sling makes it very useful for operative procedures on the lower regions of the patient's body, such as gynaecological or urology procedures.

**Material**

Polyester  
Medium . . . . . 284855  
Large . . . . . 284865  
Lifting straps incl. snap hooks. . . . 284801  
Supplied in a box of 10 slings

36



**Disposable  
Horizontal sling standard, 350 kg**

Supports the whole body.

**Suitable for**

- lifting a patient in a horizontal position.

**Material**

Polyester  
. . . . . 284631  
Supplied in a box of 10 slings

Disposable slings



37

**Disposable Leg sling II**

Used when legs need lifting and support. Possible activities are dressing, apply therapeutic sock/hose, access to hygiene, support of leg to apply ROM devices. Same design as Guldmann Leg sling.

**Material**

Polypropylene  
Supplied in a box of 10 slings . . . . .286501



38

**Disposable Multi support sling**

Supports lifting part of the body, upper body, pelvis, legs or arms.

**Suitable for**

- situations where parts of the body need to be lifted or supported – e.g. for placing X-RAY cassette.

**Material**

Polyester  
.....284223

39

**Disposable Twin Turner II**

To be used when turning a user into lateral position. Enables placing of sling under the patient without having to roll or move them. Smaller surface area of material for this sling allows access to patient's skin for inspection, treatment and care. Same design as Guldmann Twin Turner.

**Material**

Polypropylene  
Regular . . . . .287501  
Large . . . . .287511  
Supplied in a box of 10 slings



Slings and lifting accessories

Sizes and fabrics



## Sizes and fabrics

ABC Slings		Kids	Kids	Kids	Kids	Kids	XS	S	M	L	XL	2XL	3XL
		2 - 4	4 - 6	6 - 10	10 - 14	10 - 16							
<b>ACTIVE</b>													
<b>Active Micro</b>	polyester						284031	284041	284051	284061	284071	284081	
<b>Active Micro Plus</b>	polyester		281001	281011	281021		281531	281541	281551	281561	281571	281581	
<b>Active Trainer</b>	polyester		283001	283011	283021		283031	283041	283051	283061	283071	283081	
<b>BASIC</b>													
<b>Basic Basic</b>	polyester		270001	270011	270021		270031	270041	270051	270061	270071	270081	270091 <sup>x)</sup>
<b>Basic Basic</b>	net		270103	270113	270123		270133	270143	270153	270163	270173	270183	
<b>Basic Low</b>	polyester						271031	271041	271051	271061	271071	271081	
<b>Basic Low</b>	net						271133	271143	271153	271163	271173	271183	
<b>Basic High</b>	polyester		272001	272011	272021		272031	272041	272051	272061	272071	272081	
<b>Basic High</b>	net		272103	272113	272123		272133	272143	272153	272163	272173	272183	
<b>Basic Shell</b>	net						275032	275042	275052	275062	275072		
<b>Basic Comfort High</b>	polyester						277031	277041	277051	277061	277071		
<b>Basic Comfort High</b>	net						277032	277042	277052	277062	277072	277082	
<b>CUSTOM</b>													
<b>Sit-On II</b>	net						297031	297041	297051	297061	297071		
<b>Sit-On High II</b>	net						298031	298041	298051	298061	298071	298081	298091
<b>Sit-On Comfort</b>	Spacer		293001	293011			293031	293041	293051	293061	293071		
<b>Sit-On Comfort High</b>	Spacer		294001	294011			294031	294041	294051	294061	294071	294081	294091
<b>Sit-On Comfort High Hygiene</b>	Spacer						294131	294141	294151	294161	294171		
<b>Custom Amputee</b>	polyester						290031	290041	290051	290061	290071		
<b>Active Vest Kids<sup>xx)</sup></b>	polyester	2831000	283101	283111		283121							
<b>Vest for Standing Shell<sup>xx)</sup></b>	polyester			283511		283521	283531						


**Max lifting capacity for all Guldmann ABC slings: 255 kg**

x) Max lifting capacity for Basic Basic 3XL: 500 kg

xx) Max lifting capacity for Active Vest Kids / Vest for standing Shell: 100 kg

Disposable Sling		Kids	XS	S	M	L	L/XL	XL	XL/2XL	2XL	2XL/3XL	3XL
		6 - 10										
<b>Disposable Micro Plus</b> <i>Box of 5 slings</i>	polyester		281535	281545	281555	281565		281575		281585		
<b>Disposable Gait Trainer</b> <i>Box of 5 slings</i>	polyester			283545	283555	283565		283575		283585		
<b>Disposable High</b> <i>Box of 10 slings</i>	polypropylene	283615	283635	283645	283655	283665		283675		283685		
<b>Disposable High Bariatric</b> <i>Box of 10 slings</i>	polyester/cotton						283672		283682		283692	
<b>Disposable Comfort High</b> <i>Box of 5 slings</i>	polyester		277035	277045	277055	277065		277075		277085		

For max lifting capacity see the specific disposable sling description page 48-49.



Sizes and fabrics



① **Centre of sling marking**

② **Information-label**

Indicates all that the carer needs to know



③ **Label with size indication**

Makes it easy and simple to find the correct size of sling

④ **Padding in back section**

For added comfort on models with small back surface

⑤ **Positioning pocket**

All Basic slings are provided with a pocket to ease correct positioning of sling

⑥ **Coloured loops**

The coloured loops ensure symmetry in the lift

⑦ **Material**

Dirt repellent polyester characterized by its strength and elasticity

⑧ **Auxiliary straps**

It is possible to use auxiliary straps


⑨ **Padding on legs**

All slings with divided leg supports have inset foam pads on legs. Easier to put on – smoother on skin

Sizes and fabrics

Guidelines on choosing the size of a sling

All dimensions in cm



Size / cm	W1	W2 <sup>x)</sup>	W2 <sup>xx)</sup>	H
Kids 2-4			46 - 54	
Kids 4-6	25 - 30	50 - 90	50 - 65	45 - 60
Kids 6-10	28 - 33	60 - 95	51 - 72	55 - 70
Kids 10-14	31 - 36	65 - 100		68 - 80
Kids 10-16			56 - 75	
XS	33 - 37	70 - 105		75 - 82
S	36 - 40	75 - 110		80 - 87
M	39 - 43	80 - 120		85 - 92
L	42 - 45	90 - 130		90 - 97
XL	45 - 49	100 - 140		95 - 102
2XL	48 - 54	110 - 150		95 - 102
3XL	54 - 58			95 - 102
4XL	57 - 63			95 - 102
5XL	62 - 73			95 - 102

W2 = Chest measurements


x) Guideline for choosing Active Trainer and Gait Trainer

xx) Guideline for choosing Active Vest Kids

Fabrics



Polyester



Net Polyester Black



Net Polyester Grey



Spacer Polyester (Comfort)




Polyester (Disposable)



Polypropylene (Disposable)



Polyester/Cotton (Disposable bariatric)



TENCEL® (Antibacterial)

Care labelling

ABC slings



Normal washing at the indicated temperature



Do not use bleaching agent



Tumble dry at low temperature



Do not iron

Disposable sling

DO NOT WASH



All slings are provided with a label saying "DO NOT WASH" which becomes visible with a warning "DO NOT USE" when in contact with liquid as an additional safety measure. Disposable slings are used in situations where hygiene is important and to minimise the risk of cross infection.



40



**Positioning sling**

The positioning sling system makes it possible to easily turn the patient to one side or to a prone position, e.g. on the OR table. Using the sling together with the hoist prevents shear force and friction on the patient's skin and improves the working environment for the care staff.

Side positioning sling . . . . .284226

*Prone positioning package:*

*2 x slings*

*2 x Multi Support slings*

*1 x security strap*

. . . . .284225



Turning and handling options



41



**Repositioning sling, 500 kg**

A sling for use as part of the linen make-up of the patient bed; leave positioned underneath patient, ready to use for care tasks.

**Suitable for**

- rolling to side lying • repositioning patient up in bed • sitting forward head and trunk
- supine lateral transfers to another surface such as beds + stretchers • full lift in supine to change linen, or other.

**Material**

Standard polyester (photo)  
 1830 mm . . . . . 28465  
 Standard net, 1830 mm. . . . . 284658  
 Short polyester, 1350 mm . . . . . 284653  
 Long, 2130 mm . . . . . 284651  
 Supplied in a box of 10 slings



42



**Repositioning sling, 255 kg**

A sling for use as part of the linen make-up of the patient bed; leave positioned underneath patient, ready to use for care tasks.

**Suitable for**

- rolling to side lying • repositioning patient up in bed • sitting forward head and trunk
- supine lateral transfers to another surface such as beds + stretchers • full lift in supine to change linen, or other.

**Material**

Spacer, 1990 mm . . . . . 284659  
 Supplied in a box of 5 slings



43



**Repositioning sling, TENCEL®, 375 kg**

A sling for use as part of the linen make-up of the patient bed; leave positioned underneath patient, ready to use for care tasks.

**Suitable for**

- rolling to side lying, repositioning patient up in bed • sitting forward head and trunk
- supine lateral transfers to another surface such as beds + stretchers • full lift in supine to change linen, or other.

**Material**

TENCEL® is an organic material with antibacterial and temperature regulating characteristics. The special TENCEL® fibres are exceptionally soft and gentle against sensitive skin.  
 2130 mm . . . . . 284657  
 Supplied in a box of 10 slings



Turning and handling options



44

**OR Sling – Upper body**

Supports the upper body when moving to and from the operating table in a horizontal position.

**Suitable for**

- the length and design of the sling makes it very useful for operative procedures on the lower regions of the patient's body, such as gynaecological or urology procedures.

**Material**

Polyester

Medium . . . . .	284851
Large . . . . .	284861
Lifting straps incl. snap hooks. . . . .	284801

45

**Multi support sling**

Supports lifting part of the body, upper body, pelvis, legs or arms.

**Suitable for**

- situations where parts of the body need to be lifted or supported – e.g. for placing X-RAY cassette.

. . . . . 28467



46



**Turner**

Enables placing of sling under patient without having to roll or move them. Smaller surface area of material for this sling allows access to patient's skin for inspection, treatment and care.

. . . . . 28700

47



**Twin Turner**

Makes it possible to roll a patient on both sides without removing the sling. Enables placing of sling under patients without having to roll or move them. Smaller surface area of material for this sling allows access to patient's skin for inspection, treatment and care.

Regular . . . . .	28750
Large . . . . .	28751

48



**Leg sling**

Used when legs need lifting and support. Possible activities are dressing, apply therapeutic sock/hose, access to hygiene, support of leg to apply ROM devices.

. . . . . 28650  
*Supplied in a box of 10 slings*

Lifting hanger options



49



**Lifting hanger, medium, 375 kg**

Suitable for persons using a lifting sling in sizes M-3XL.

**Function**

To be used in connection with uncomplicated lifts and when using standard slings.  
 .....556890  
 With adjustable resistance turning .556930

50



**Lifting hanger, small, 375 kg**

Suitable for persons using a lifting sling in sizes XS and S.

**Function**

To be used in connection with uncomplicated lifts and when using standard slings.  
 The lifting hanger is recommended when using lifting slings size S and down.  
 .....556880  
 With adjustable resistance turning .556920

51



**Lifting hanger, x-small, 375 kg**

Suitable for persons using a lifting sling in XS and Kid Sizes.

**Function**

To be used in connection with uncomplicated lifts and when using standard slings.  
 The lifting hanger is recommended when using the "kids" lifting slings.  
 .....556870  
 With adjustable resistance turning .556910

52



**Cross hanger, 500 kg**

Suitable when larger spacing between back and leg straps of the lifting sling is required. To be used when lifting with GH3+ and GH3 Twin ceiling hoist.

The attachment top comes in two different versions:  
 GH3+ 1 strap ..... 550444  
 GH3 Twin 2 straps ..... 550800


53



**Cross hanger, 400 kg**

Suitable when larger spacing between back and leg straps of the lifting sling is required. To be used when lifting with GH1, GH3, GH3+, GH3 Twin ceiling hoist and GL5 mobile lifter.  
 .....561610

Lifting hanger options



54

**H Hanger, 375 kg**

Suitable when lifting persons where a large spacing between back and leg straps of the lifting sling is required, e.g. when lifting persons with a large abdomen/ stomach. The H Hanger is available with side hangers in 3 different sizes.  
 H Hanger – centre hanger. . . . . 556950

**Note:** In addition you need to order 2 hangers of either x-small (556870) small (556880) or standard (556890)

55



**Swing kit**  
 ..... 554903

56



**Connecting bar**  
 Optional for GH3 Twin so diverse hangers can be utilized  
 ..... 550544



57



**Snap hook**  
 ..... 561628

Slings and lifting accessories

Horizontal lifters



Horizontal lifters



58

**Horizontal lifter, foldable, 255 kg**

Suitable when a horizontal and stable lift is required, with weight adjustment. Can be used with either a complete lifting sheet or straps.

Horizontal lifter

incl. one lifting sheet . . . . . 28445

Lifting sheet . . . . . 284485

Set of lifting straps incl. stiffener . . . . 28448

Extra clamp . . . . . 2844529

59

**Horizontal lifter, foldable, 255 kg  
2 strap version**


Suitable when a horizontal and stable lift is required, with weight adjustment. To be used with GH3 Twin. Can be used with either a complete lifting sheet or straps.

Horizontal lifter

incl. one lifting sheet . . . . . 28444

Lifting sheet . . . . . 284485

Set of lifting straps incl. stiffener . . . . 28448



60

**Trolley for horizontal lifter**

. . . . . 28443



Horizontal lifters

61

**Horizontal lifting support, 375 kg incl. swivel link**

Horizontal lifting support with infinitely variable weight adjustment. Can be used with GH3 ceiling hoists or with Guldmann mobile lifter.

- ..... 28456
- Horizontal sling, standard ..... 28463



62



**Horizontal lifting support, 375 kg**

Horizontal lifting support with 10 different lifting positions. Can be used with GH3 ceiling hoists or with Guldmann mobile lifter.

- ..... 28466
- Horizontal sling, standard ..... 28463

63



**Mortuary lifting straps**

Ease lifting in mortuary. Can take larger corpses / bodies. The Individual straps allow easy access.

- Standard incl. stiffener ..... 284640

64



**Bath chair – rustproof, 255 kg**

Suitable when a stable lift into a swimming pool is required.

- ..... 14600

65



**Bath stretcher – rustproof, 255 kg**

Suitable when a stable and horizontal lift into a swimming pool is required.

- ..... 14610

Accessories



66

**Head support**

Supports the head without stimulating the neck reflexes. Suitable for persons with no head control.

.....28325

67



**Sheepskin insets**

Suitable for persons who need more gentle support of the thighbone

.....28315

68




**Cover for lifting hanger**

XS .....28316

S .....28317

M .....28318

69



**Marking clips, 4 pcs**

Red .....28477

Blue .....28478

Green.....28479

70



**Auxiliary straps**

Ease correct positioning of the user in the wheelchair. Use auxiliary straps to attach for handling.

.....28355

71



**Extension straps for slings**

Used in situations to achieve various sitting positions or where different lifting heights are required. *Set of two pcs.*

20 cm. ....28470

30 cm. ....28472

72



**Lifting strap extension**

0.5 m .....550459-500

1.0 m .....550459-1000

1.5 m .....550459-1500

2.0 m .....550459-2000



Accessories



73



**Smart Glide – rolls of 100 m**


A slide sheet used for reducing friction between the user and the surface.  
Maximum user weight 255 kg.

**Material**

Nylon with silicone coating.

Width M 100 cm. . . . .284472  
Width L 140 cm. . . . .284473

74



**Smart Glide Cutter**

. . . . .284479

75



**Smart Glide Dispenser**

Smart Glide Dispenser for wall mounting in either vertical or horizontal position.  
Flexible for different width of rollers.

. . . . .561470

76

**Smart Glide Precut**

Each piece of Smart Glide Precut is in a separate bag which can be used for storage. The box is prepared for wall mounting. Maximum user weight 255 kg.

100 cm x 125 cm  
5 boxes of 20 pcs. . . . .561500  
1 box of 20 pcs . . . . . 561500-1

100 cm x 145 cm  
5 boxes of 20 pcs. . . . .561501  
1 box of 20 pcs. . . . . 561501-1



Accessories

77

**Digital scale for ceiling hoist or mobile lifter**

Battery operated digital scale.

- Digital scale, 300 kg
- Standard Model RS . . . . .560602
- Class III Model RS, verifiable . . . . .560604
- Quick connect . . . . .560648



78



**Adapter kit – horizontal lifting devices**

. . . . .560633

79



**Swivel link**

For horizontal lifters and horizontal lifting supports.

. . . . .284520

80



**Holder for horizontal lifting support**

. . . . .556968



81



**HC strap rest**

To be placed on the lifting strap of the hoist

- 1 pcs . . . . .559516
- 10 pcs . . . . .559943

82



**Holder for lifting hanger**

. . . . .556940

Accessories



83



**Training kit for Positioning lock**

For the Guldmann Positioning lock a complete training kit is available containing various training accessories

..... 293560


The kit contains:

- A) Box with 30 m training tube: Orange (light) ..... 558136
- B) Box with 30 m training tube: Lime green (medium) ... 558137
- C) Box with 30 m training tube: Plum (hard) ..... 558138
- D) 2 hand grips for training tube ..... 558139
- E) Snap hook ..... 561628
- F) 2 connecting straps ..... 29356
- G) 2 elastic lines: Blue (light) and black (hard) ..... 561629
- H) Multi support sling ..... 28467

A-H can be ordered separately



| Time to care |



GB

**V. Guldmann A/S**  
Graham Bells Vej 21-23A  
DK-8200 Aarhus N  
Tel. +45 8741 3100  
info@guldmann.com  
www.guldmann.com

**Australian Distributor**  
HLS Healthcare Pty Ltd  
67 Henderson Road,  
Rowville, Victoria, 3178  
Tel. 1300 931 893  
sales@hlshealthcare.com.au  
www.hlshealthcare.com.au



© Guldmann AUS-2799\_2/10/2023

**Guldmann™**